

Holistic trichology is the awareness of having applied all the best

3Active raises the effectiveness of integrators to a power, driving amino acids, minerals, nourishing and curing substances right into the bulb.

Massage, ungluing, phototherapy multiply their effectiveness while simultaneously applied by 3Active. Active principles and shampoo complete the holistic treatment.


3 active with its triple action is the basis of the holistic trichology. Photostimulation for the revitalization of the bulb.

Vibrational massage for the superficial circulation.  
Ungluing for the deep circulation.



Photostimulation for the revitalization of the bulb.

Vibrational massage for the superficial circulation.



Ungluing for the deep circulation.

# 3Active works



3 scientific solutions in a single device.

[www.3Active.it](http://www.3Active.it)

# 3active

3active acts with 3 simultaneous actions, which raise its effectiveness to the third power.

## 1 Deep massage



3Active drives the blood towards the upper part of the skull.

The deep massage creates the necessary space to let the blood get into each bulb of the hair. In doing so, it nourishes it and it lets the exhausted and intoxicated blood, which would poison the hair, flow outside - similarly to a sponge squeezed and washed under Running water.



Fresh blood, full Of active principles, regenerates the life of the hair from its inside. 3active conveys its highest effectiveness when used with active principles or integrators, initiating a holistic approach for the care and regrowth of the hair.


## 2 VIBRATIONAL MASSAGE



The activation of the circulation is amplified by a vibration with variable frequency.

A program starts the activation when 3active is pushing the blood towards the top of the head, removing exhausted blood and toxins stored because of poor circulation.

## 3 PHOTOSTIMULATION



Low level light therapy. With this search key you can find 7 million pages. It is the most popular technology in this moment. 3active contains fifteen 633 nm sources. 3active multiplies its effects photostimulating the scalp preactivated by a deep and superficial massage action.

## ACTIVE PRINCIPLES



With its patented deep massage, 3active is unbeatable in helping the deep penetration of active principles.

### Only applied scientific researches.

The benefits of the massage called 卍 gungling卍 h are already known, the photostimulation is a new technique. We have invented a series of technical solutions to apply simultaneously all therapies which are recognized as valid by common experience and scientific research. 3 instruments in a single device.

*In the following list, just some of the scientific researches which guided us in the realization of 3active.*

### **Researches about photostimulation:**

- 1) *Combined 633-nm and 830-nm LED treatment of photoaging skin. From: Journal of Drugs in*
- 2) *Effect of pulsing in low-level light therapy Hashmi JT, Huang YY, Sharma SK, Kurup DB, De Taboada L, Carroll JD, Hamblin MR.*

*1)Journal of Investigative Dermatology, Vol 92, 725-726,*

*2) PMID: 8628793 [PubMed - indexed for MEDLINE]*

*P.Klemp, K.Peters and B.Hansted Dipartimento di dermatologia, Università di Copenhagen,*